


# E-TRAVELLER PURSUIT

PORTABLE MOBILITY SCOOTER • PS3500



**SOLID PUNCTURE  
PROOF TYRES**



**FOLDS FOR EASY  
TRANSPORT**



**FRONT AND REAR  
SUSPENSION**



# E-TRAVELLER PURSUIT

## PORTABLE MOBILITY SCOOTER • PS3500

Back to due popular demand! The Pursuit carries the ultimate in features for any travel enthusiast. Its lightweight design and tight turning circle make it ideal for taller individuals, allowing it to remain agile and easily manoeuvrable. Each Pursuit has a drink bottle holder and two front basket options, that can be carried together or individually. Choose from a collapsible front basket or an insulated front pouch, or opt to place the front pouch into the front basket for the best out of two worlds. With solid puncture-proof tyres and bright LED headlights, you'll never have to worry about getting a flat or losing your way in the dark. With up to 15km of range per charge and speed of up to 6km/h you can easily travel long distances in a short amount of time.

### SPECIFICATIONS

<b>Dimensions</b>	940x470x850-910mm	<b>Ground Clearance</b>	38mm
<b>Total Weight</b>	27.1kg	<b>Turning Radius</b>	1100mm
<b>Heaviest Part</b>	20kg	<b>Slope Ability</b>	10°
<b>Drive Wheel</b>	178mm	<b>Battery Range</b>	Up to 12km
<b>Motor</b>	120W	<b>Battery</b>	24V 6Ah Lithium
<b>Weight Capacity</b>	135kg	<b>Seat Size</b>	420x380mm
<b>Tyres</b>	Solid - Puncture Proof	<b>Anti-Tip Mechanism</b>	Rear Anti-Tip Wheels