

E-TRAVELLER VOYAGER

REMOTE CONTROL WHEELCHAIR • EW4330



**EASY TO USE
JOYSTICK CONTROL**



**FOLDS FOR EASY
TRANSPORT**



**LIGHTWEIGHT
ALUMINIUM FRAME**



E-TRAVELLER VOYAGER

REMOTE CONTROL ELECTRIC WHEELCHAIR

The E-Traveller Voyager is a lightweight folding electric wheelchair, which is an ideal companion for everyday travel. This compact wheelchair can be folded within seconds for easy storage and transportation in planes, cars, trains and more. It is equipped with an LCD digital panel, remote control operation, puncture-proof tyres and a lightweight aluminium frame.

SPECIFICATIONS

Weight Capacity	100kg	Seat Dimensions WxD (mm)	460x420
Dimensions LxWxH (mm)	1040x620x950	Max Slope	12°
Folded Dimensions (mm)	340x610x750	Tyres	Solid Puncture-Proof
Max Speed	6km/h	Drive Wheel	12 inch
Battery Range	12km	Motor	180W x 2
Battery	24V 6Ah	Control	Joystick and Remote
Weight Without Battery	17.5kg	Total Weight	21kg



**SCAN OR CLICK
TO FIND OUT MORE**