## HEARTWAY MIRAGE LUX

MEDIUM SIZED MOBILITY SCOOTER • MS0020









UPGRADED LUXURY COMFORT SEAT



REAR SUSPENSION





## **HEARTWAY** MIRAGE LUX

## **MEDIUM SIZED MOBILITY SCOOTER**

The Heartway Mirage Lux sees the addition of a luxury comfort seat to go along with all the features seen on the Mirage. Perfect for running errands at the local shops and exploring the great outdoors alike. It boasts a travel range of up to 30km per charge, large chunky tyres and rear suspension to keep you comfortable while you out and about. The display panel has been designed with ease of use in mind, featuring large soft touch buttons for intuitive operation. The seat can be adjusted in height, allowing for a 360-degree swivel for easy transfer.

SPECIFICATIONS			
Weight With Batteries	87kg	Tyres	305mm
Suspension	Rear Suspension	Dimensions LxW (mm)	1400 x 620 x 1280
Weight Capacity	135kg	Swivel Seat	Yes
Adjustable Tiller Angle	Yes	Adjustable Headrest	Yes
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Adjustable Armrests	Yes
Battery	12V/50ah x 2	Drive Wheel	305mm
Battery Range	Up To 30km	Ground Clearance	130mm
Maximum Speed	10km/h	Adjustable Seat Height	600mm-680mm
Motor	550W 4-Pole	Lights	Yes
High/Low Speed Selection	Yes	Turning Circle	1062mm

