

QUINGO CONNECT

ELECTRIC WHEELCHAIR



**SOLID PUNCTURE
PROOF TYRES**



**EASY TO USE
JOYSTICK CONTROL**



**TWO STORAGE
BAGS**



QUINGO CONNECT

ELECTRIC WHEELCHAIR

The Quingo Connect has a fantastic folding design and uses a simple push-pull folding frame. In under 5 seconds, the Quingo Connect can be transformed into a compact size ideal for easy storage, for transportation or can be pulled along. The Quingo Connect is extremely maneuverable and easy to control via its simple to use joystick control. The Connect can turn on its own axis making it easy to navigate tight and confined spaces.

SPECIFICATIONS

Weight With Batteries	27.3kg	Tyres	Solid - Puncture Proof
Battery Weight	2.8kg	Dimensions LxWxH (mm)	970x590x950
Heaviest Part	24.5kg	Folded Dimensions (mm)	590x370x810
Weight Capacity	114kg	Distance Between Armrests	470mm
Controller	Joystick	Drive Wheel	178mm
Cushion	Contoured	Ground Clearance	175mm
Battery	24V 11.4ah	Seat Dimensions	460mm x 420mm
Battery Range	Up To 15km	Seat Height From Ground	520mm
Maximum Speed	6.4km/h	Slope Ability	3°
Motor	70W	Turning Circle	825mm
Number of Components	1 + Battery	Storage	Back pocket & under seat bag



www.wholesalemobility.com.au