

QUINGO PLUS

5 WHEEL MOBILITY SCOOTER



**INNOVATIVE
5-WHEEL DESIGN**



**LARGE
COMFORTABLE SEAT**



**KERB CLIMBING
TECHNOLOGY**



QUINGO PLUS

5 WHEEL MOBILITY SCOOTER

The Quingo Plus is an upgraded model of the Quingo classic, featuring everything we love about the Classic and the 5-wheel Quingo design plus additional comfort and usability features, such as high and low-speed selection, a broader range and a large adjustable captain seat. The Plus is classed as a medium-sized mobility scooter but can outperform many other large gophers. However, thanks to the 5-wheel design, it still allows for the same maneuverability that you are used to from smaller 4-wheel models.

SPECIFICATIONS

Weight With Batteries	103kg	Tyres	Pneumatic
Battery Weight	-	Dimensions LxW (mm)	1270x650
Heaviest Part	-	Swivel Seat	Yes
Weight Capacity	159kg	Distance Between Armrests	450mm-630mm
Adjustable Tiller Angle	Yes	Drive Wheel	250mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	55mm
Battery	12V 50Ah x2	Seat Dimensions	480mm x 440mm
Battery Range	Up To 48km	Seat Height From Ground	600mm - 700mm
Maximum Speed	6.4km/h	Slope Ability	12°
Motor	400W 4-Pole	Turning Circle	1070mm
Number of Components	1	Suspension	Yes



www.wholesalemobility.com.au