QUINGO TOURA 2

5 WHEEL MOBILITY SCOOTER





INNOVATIVE

5-WHEEL DESIGN





DIGITAL DASHBOARD

FRONT BASKET ATTACHMENT





QUINGO TOURA 2

5 WHEEL MOBILITY SCOOTER

The Quingo Toura 2 is the top of the line scooter of the 5-wheel Quingo range. It features countless comfort and safety features, all while maintaining excellent maneuverability for such a substantial machine. With a carrying capacity of 220kg and a travel range of up to 88km this scooter makes the perfect car replacement. Featuring an easy to read digital dash with a battery meter gauge, easy to use speed dial and push-button lighting control the dash has everything that you need exactly where you need it.

SPECIFICATIONS			
Weight With Batteries	163.9kg	Tyres	Pneumatic
Battery Weight	-	Dimensions LxW (mm)	1420x700
Heaviest Part	-	Swivel Seat	Yes
Weight Capacity	220kg	Distance Between Armrests	540mm-720mm
Adjustable Tiller Angle	Yes	Drive Wheel	305mm
Anti-Tip Mechanism	No	Ground Clearance	100mm
Battery	12V 100Ah x2	Seat Dimensions	540mm x 470mm
Battery Range	Up To 88km	Seat Height From Ground	720mm - 790mm
Maximum Speed	10km/h	Slope Ability	12°
Motor	650W 4-Pole	Turning Circle	1320mm
Number of Components	1	Suspension	Yes

