QUINGO ULTRA

PORTABLE 5 WHEEL MOBILITY SCOOTER





INNOVATIVE 5-WHEEL DESIGN



USB PORT FOR CHARGING PHONE



PULLS APART
TO FIT IN CAR BOOT





QUINGO ULTRA

PORTABLE 5 WHEEL MOBILITY SCOOTER

The Quingo Ultra is a portable 5 wheel mobility scooter and the newest model in the Quingo range. The Ultra has all the trademark features of the quingo range including the 5 wheel quintell design, rear anti-tip wheels, a front LED headlight, but also has the added benefit of pulling apart to fit in your boot for transportation. One added feature that has been added to the Ultra is the handy USB phone charging port that is located on the tiller.

SPECIFICATIONS			
Weight With Batteries	-	Tyres	Pneumatic
Battery Weight	4.5kg x 4	Dimensions LxW (mm)	1050x580
Heaviest Part	17.8kg	Swivel Seat	Yes
Weight Capacity	140kg	USB Phone Charging	Yes
Adjustable Tiller Angle	Yes	Drive Wheel	220mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	50mm
Battery	4 x 15 Ah SLA	Seat Dimensions	580mm x 370mm
Battery Range	Up To 29km	Seat Height From Ground	540 - 590mm
Maximum Speed	6.4km/h	Slope Ability	6°
Motor	250W	Distance Between Armrests	460mm
Number of Components	4 + Basket	Storage	Front Basket

