CUINGO VITESS 2 5 WHEEL MOBILITY SCOOTER







GUINGO VITESS 2 5 WHEEL MOBILITY SCOOTER

The Quingo Vitess 2 is the perfect mobility solution for more hilly terrain. It features a Digital Dash, superior lighting & more power. The unique 5 wheel stability system on all Quingo's provides the excellent posture and driving position of a 3 wheeler with the stability and safer cornering dynamics of a 4 wheeler. Along with the best possible anatomical driving position all Quingo's also benefit from a staggering 80% more foot area than in a comparable length four-wheel vehicle.

SPECIFICATIONS			
Weight With Batteries	116kg	Tyres	Pneumatic
Battery Weight	-	Dimensions LxW (mm)	1420x700
Heaviest Part	-	Swivel Seat	Yes
Weight Capacity	182kg	Distance Between Armrests	450mm-630mm
Adjustable Tiller Angle	Yes	Drive Wheel	178mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	75mm
Battery	12V 50Ah x2	Seat Dimensions	470mm x 430mm
Battery Range	Up To 48km	Seat Height From Ground	640mm - 740mm
Maximum Speed	10km/h	Slope Ability	12°
Motor	400W 4-Pole	Turning Circle	1130mm
Number of Components	1	Suspension	Yes



www.wholesalemobility.com.au