

SWEETRICH AIR

PORTABLE MOBILITY SCOOTER • PS6200



**LITHIUM BATTERY &
USB PHONE CHARGING**



**PUNCTURE PROOF
TYRES**



**LIGHTWEIGHT
COMPACT DESIGN**



SWEETRICH AIR

PORTABLE MOBILITY SCOOTER

The Sweetrich AIR is a lightweight portable mobility scooter, with the heaviest part weighing just 15KG. Airline-approved with a detachable lithium battery, it ensures seamless travel. Safety is paramount with LED headlights and rear brake light. The compact front wheels allow for an incredible turning circle. Explore cities or terrains comfortably with this innovative solution that redefines mobility equipment.

SPECIFICATIONS

Weight Capacity	136kg	Total Weight	26kg
Dimensions LxWxH	1000mm x 525mm x 920mm	Max Slope	8°
Folded Dimensions LxWxH	1000mm x 525mm x 260mm	Tyres	Solid Puncture-Proof
Max Speed	6kph	Frame	Aluminium
Battery Range	Up to 15km	Motor	250W
Battery	24V 12Ah Lithium	Turning Circle	On the spot
Battery Weight	3.4kg	Suspension	Rear